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Carbon fibre produces heat thus no magnetic radiation.

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Using low electrical voltage insures much safer in using our Products.

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Yu Yang Bianstone Pad is a good helper for all drivers.

It can improve the backbone of the driver.

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The voltage of the Pad is 12V and can be used in most cars. The warmth of the Pad can remove the tiredness and sore of the back and the hip.

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Body temperature affects your health and blood circulation is influenced by the coldness and warmth of the body.

砭石熱墊可幫助人體活血益氣、疏通經絡、扶正祛邪、平衡陰陽。Bianstone Heating Pad can improve the blood circulation, enhance the body energy, eliminate toxic substances out of the body and balance Yin and Yang of the body.

遠紅外線中波長 8 ~ 14 微米，是肉眼看不到的太陽光，又稱生命之光，可滲透到身體內，提升體溫活化水分及新陳代謝，分解代謝廢物。
Bianstone possesses infra-red energy of wave length between 8 to 14 micrometer which is also called solar light or life light. This wave length of infra-red energy can penetrate into our body to enhance our body temperature and betterment our metabolism. As a result our immune system can be improved to protect our body.

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Heat Treatment: improve our cell metabolism, speed our blood movement and reduce stagnant plague in our blood vessels as well as eliminate coldness and wet.

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Infra-red wave: an amazing energy field of far infra-red of wave length between 8 to 14 micrometer can benefit our circulation of meridians, increase body tissue energy and improve our immune system.

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Bianstone possesses thirty or more rare elements such as Strontium, Titanium, Chromium, Manganese and Zinc. When these elements contact our skin, they can activate the blood circulation inside our blood capillaries.

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1. Warming the back and acupuncture points such as BL13, BL15, GV9 can relieve and ease the muscle of the back, reduce fatigue thus enhance the mental spirit of the driver.

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緩解腰部肌肉疲勞、強腰固腎，對長期駕駛者的腰骶酸痛，有明顯效果。
2. Warming the waist and acupuncture point GV4 can strengthen the muscle of the waist and increase Yang Qi (solar energy). For drivers who suffer chronic back pain can feel a significant comfortable effect.

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3. Warming the hip can ease and relieve the stress muscle and tendon of the hip. For long distance driver, Yu Yang Warm Pad is a very good companion

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Place the plug into the cigarette switch after starting the car. Press the \odot for a few seconds to start the operation. After switching on the Pad, wait for 10 minutes until the default temperature reach 50℃ with a difference +5℃ or -5℃. The default fixed time is 60 min. and the count-down begins. Press the \odot for a few seconds to turn off the Pad. Remove the plug from the electric power when not in use.

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Count-down setting: Press the \odot to adjust the time duration. There are three count-down time displays which are 60/90/120 minutes. Press the “ \oplus ” switch or “ \ominus ” to adjust the time required.

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Temperature setting: Press the \odot to adjust the temperature. There are 5 temperature settings which are 40/45/50/55/60℃. Press the “ \oplus ” switch or “ \ominus ” to adjust the temperature required.

4. 聲音提示：倒數計時完成，會有一聲“滴”，有故障一直長“滴”。
Sound indicator: 1 beep sound means the count-down reaches zero. Once the Heating Pad is not functioning well, a long beep sound will be heard.

5. 液晶顯示 E1、E2、E3 表示有故障，請聯絡本公司維修。
When the display panel shows E1, E2 or E3, it means the Heating Pad needs repairing. Please contact Our Company immediately.

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Precautions

- 心臟病或心律不正者慎用。
Use with great care for those people who are suffering from heart disease or irregular heart pulses.
- 女性生理期間或懷孕請勿在腹部使用。
Do not apply on abdomens during menstruation period or pregnancy.
- 酒後禁止使用。
Do not use the equipment after drinking wine.
- 若本產品跌落或遭撞擊而受損時，請停止使用。
Stop using the equipment if it is damaged.
- 本機若出現故障，請交本公司修理，若由非專業人士修理，可能因操作不當而導致意外等。
In case of malfunction of the product, please take it back to our Company for repairing. Never attempt to repair it by non-professional persons to prevent unexpected accidents.
- 請勿將本產品在接近水源的地方使用，例如浴室、浴缸、臉盆或其它高度潮濕的地方，會引致故障。
Never place the equipment near wet area such as bathroom, basin or other high humidity positions which may cause damage to the equipment.
- 汽車熄火後，請將電源插頭拔出，以防走電。
After the switch of the car is turned off, please immediately pull out the plug to avoid leakage of electricity.

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Cross Section Description